



ALIGNING SIGNALS

~ configuring the data ~




- ② shift BM/MC signal BACK by ~5 hours
- account for time zone differences
 - account for machine drift differences




- ③ put NIRS and BM/MC signals onto the same time zone
- NIRS sampling rate = 15 points / minute \therefore every 4 seconds
 - BM/MC sampling rate = 60 points / minute \therefore every 1 second

FINDING HYPOXIC EPISODES

~ extracting useful information from data ~



- ④ Get the start of each hypoxic episode by identifying them in the SpO₂ signal
- to do this, though, I filtered SpO₂ values to retrieve all those $\leq 80\%$



this is done by assigning a 1 to a low SpO₂ value that was preceded by a low value

a break in a sequence of 1s indicates the end of an episode

distinguish

and you can count the # of 1s to determine episode length
(because the sampling rate of the SpO₂ signal is 1 point/second)

⑤ only keep episode lengths of 20 seconds or longer

filter

and store the start of each hypoxic episode & its duration (in seconds)

GET SpO₂ & rSO₂ DATA DURING HYPOXIAS

~ matching the datasets ~

⑥ the signals generated for the rSO₂ data is at a lower sampling rate, so obtaining a match in time intervals will be tricky

a_1

a_2

a_3

a_4

BM/MC

NIRS

for each of the hypoxia start times (a_1, a_2, a_3, a_4) get the absolute value difference in times between it and every single recorded time in the NIRS signal

EXAMPLE

matrix form

$[b_1 \ b_2 \ b_3 \ b_4 \ b_5 \ b_6 \ \dots \ b_{i-1} \ b_i]$
time points w/ difference of 4 seconds between consecutive points

$[b_1 \ b_2 \ b_3 \ b_4 \ b_5 \ b_6 \ \dots \ b_{i-1} \ b_i]$
- $[a_1 \ a_1 \ a_1 \ a_1 \ a_1 \ a_1 \ \dots \ a_1 \ a_1]$

↓ is essentially subtracting a row matrix by a scalar to produce a row matrix of durations


$[|b_1 - a_1| \ |b_2 - a_1| \ |b_3 - a_1| \ \dots]$

and get the smallest value in this matrix

if this value $|b_k - a_i| \leq 4$ seconds, then we've found a match

→ if this value were > 4 seconds, that means we don't have NIRS data captured during the hypoxic episode, so we'll "discard" it

⑦ we'll keep ± 1 minute before the start and after the end of the hypoxia for context



⑧ plot each of the hypoxia intervals